



# Versana Essential™

## Transducer Guide



	Applications	Footprint	Frequency Range	FOV	Biopsy Guide
 4C-RS	Abdominal, OB,GYN, Vascular, Urology, Thoracic, Pediatric, MSK, Interventional Guidance	18.3 x 66.2 mm	2.0 – 5.0 MHz	58°	E8385NA
 L6-12-RS	Vascular, Pediatric, Small Parts, MSK, Thoracic, Interventional Guidance	11.4 x 47.0 mm	4.0 – 13.0 MHz	38.4 mm	H40432LC
 E8C-RS	OB, GYN, Urology, Transvaginal, Transrectal, Interventional Guidance	16.9 x 21.2 mm	4.2 – 10.0 MHz	128°	E8385MJ or H40412LN
 E8Cs-RS	OB, GYN, Urology, Transvaginal, Transrectal, Interventional Guidance	18.6 x 13.9 mm	4.0 – 10.0 MHz	168°	E8385MJ or H40412LN
 3Sc-RS	Cardiac, Abdominal, Vascular, Transcranial, Thoracic, Interventional Guidance	23.7 x 18.4 mm	1.7 – 4.0 MHz	120°	H46222LC
 6S-RS	Cardiac Pediatric, Vascular, Pediatric, Transcranial, Interventional Guidance	23.5 X 16.8 mm	2.5 – 7.0 MHz	120°	Not available
 RAB2-6-RS	Abdominal, OB, GYN, Urology, Interventional Guidance	62.2 x 34.0 mm	2.0 – 6.0 MHz	70° (B), 70° x 84° (Volume scan)	H48681ML
 8C-RS	Pediatric, MSK, Cardiac Pediatric, Transcranial, Interventional Guidance	12.0 x 22.0 mm	4.2 – 10.0 MHz	131°	Not available
 LK760-RS	MSK, Interventional Guidance	67 x 13 mm	3.5 – 10.0 MHz	59.9 mm	Not available

© 2022 General Electric Company.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Versana Essential are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

November 2022  
JB22814XX  
DOC2799155

